



918 E. Divide Avenue
Bismarck, North Dakota 58501
Phone: (701) 328 -5150
Fax: (701) 328-5200

AUTHORIZATIONS

Title	Name	Signature
SOP Author	McKenzie Schick	
Program Manager	Aaron Larsen	

QUALITY CONTROL/QUALITY ASSURANCE DOCUMENTATION

Title: Measuring Dissolved Oxygen by the “Winkler-Azide” Method
 Type: Standard Operation Procedure #7.01
 Version: 3.0
 Date: 01/20/2021
 Author: McKenzie Schick

REVISION HISTORY

Revision	Change Description	Date	Authorization
1	Reviewed by AL	11/16/21	AL

ACKNOWLEDGEMENTS

(Place to acknowledge peer reviewer)

THIS PAGE LEFT INTENTIONALLY BLANK

TABLE OF CONTENTS

1.0 SCOPE AND APPLICABILITY	2
2.0 SUMMARY OF METHOD	2
3.0 HEALTH AND SAFETY WARNING	2
4.0 CAUTIONS.....	2
5.0 INTERFERENCES	2
6.0 FIELD PERSONNEL QUALIFICAITONS/RESPONSIBILITIES	3
7.0 EQUIPMENT AND SUPPLIES	3
8.0 PROCEDURE	3
9.0 DATA AND RECORD MANAGEMENT	4
10.0 QUALITY ASSURANCE AND QUALITY CONTROL	4
11.00 REFERENCES.....	4
APPENDIX A. Field Reporting Form	
APPENDIX B. SOP Acknowledgement and Training Form	

1.0 SCOPE AND APPLICABILITY

This document presents the North Dakota Department of Environmental Quality, Division of Water Quality's (DWQ) Standard Operating Procedure (SOP) for measuring dissolved oxygen (DO) by the "Winkler-Azide" method. This SOP applies to all DWQ field staff, non-DWQ cooperators, and citizen volunteers.

2.0 SUMMARY OF METHODS

Measurements made using the "Winkler-Azide" method are considered the most accurate measurement of dissolved oxygen. The "Winkler-Azide" method can be used for making direct measurements of dissolved oxygen in surface waters, field calibrating a dissolved oxygen meter, and performing regularly scheduled dissolved oxygen meter calibrations. The following are the procedures for measuring dissolved oxygen by the "Winkler-Azide" method. A detailed description of the method can be found in Standard Methods, (17th Edition, APHA, 1989).

3.0 HEALTH AND SAFETY WARNING

In field situations, personnel should take appropriate precautions when operating watercraft and working on, in, or around water. All boats used for sampling should be equipped with safety equipment such as personal flotation devices (PFD's), oars, air horn, etc. North Dakota's boating laws and rules shall be followed by all field personnel.

In the laboratory setting, appropriate attire should be donned such as protective eyewear, gloves, and aprons in order to minimize any potential splatter associated with this lab process.

4.0 CAUTIONS

Care should be taken when collecting the sample and when filling the BOD bottle so as to not introduce air bubbles into the sample.

5.0 INTERFERENCES

The most common interferences to the "Winkler-Azide" method are oxidizing and reducing agents, nitrate ion, ferrous iron, and organic matter.

6.0 PERSONNEL QUALIFICATIONS/RESPONSIBILITIES

All personnel taking field measurements must read this SOP annually and acknowledge they have done so via a signature page (see Appendix B). New field personnel must also demonstrate successful performance of the method. The signature page will be signed by both trainee and trainer to confirm that training was successfully completed and that new personnel is competent in carrying out this SOP. The signature page will be kept on-file at DWQ along with the official hard copy of this SOP.

7.0 EQUIPMENT AND SUPPLIES

- 300 ml BOD bottles
- 300 ml manganous sulfate solution
- 300 ml alkaline-iodide-azide reagent
- Concentrated sulfuric acid
- Graduated cylinder capable of measuring 203 ml
- Titration equipment capable of measuring to 0.1 ml
- 0.025 M sodium thiosulfate solution
- Starch indicator solution
- Kemmerer or Van Dorn type sampler

8.0 PROCEDURE

1. Obtain a sample by completely filling a 300 ml BOD bottle. The bottle is filled slowly taking care to exclude air bubbles. Note: When measuring dissolved oxygen for calibration purposes two samples should be collected at the same time and place. The results of the two should agree within +/- 0.2 mg/l. If not, repeat the procedure.
2. Insert the glass stopper to remove excess water trapped in the neck of the bottle.
3. Add 2 ml of manganous sulfate.
4. Add 2 ml of alkaline iodide-azide
5. Stopper the bottle, taking care to exclude air bubbles.
6. Mix well by inverting and righting the bottle rapidly for 1 minute, at which time a flocculent will have formed.
7. Let the flocculent settle until approximately 1/3 of the bottle has cleared, then mix again.
8. Let the flocculent settle until approximately 2/3 of the bottle has cleared.

9. Add 2 ml sulfuric acid, stopper, and mix. Note: Take care when adding sulfuric acid as any spillage can cause severe burns. If any acid is spilled it should be flushed immediately with water.
10. At this point the oxygen within the sample is fixed and the sample may be stored in the dark for up to 8 hours.
11. Using the graduated cylinder, transfer 203 ml of the sample to the titration vessel.
12. While gently swirling sample, titrate with the sodium thiosulfate solution to a pale straw color.
13. Add 2 ml starch indicator solution. Sample will turn blue.
14. Continue to titrate and swirl the sample until the first disappearance of blue color.
15. Record to nearest 0.1 ml the volume of titrant used as mg/l dissolved oxygen in the sample.

9.0 DATA AND RECORDS MANAGEMENT

Sample information will be recorded on the field form (Appendix A). Once personnel reach the office, data recorded on the field form are entered into the DWQ Sample Identification Database (SID). Field notes should be used to record any quality control activity performed such as measurements taken by more than one sampler, or to record any sampling conditions that may have interfered with the reading such as high winds/wave action, cattle in water, observed flow, water surface, water clarity, water color, water odor, visual algae cover, number of dead fish, present weather, estimated inches of rain fall in past 72 hours, and any comments. Field forms and notes should be stored in the appropriate project folder at DWQ.

10.0 QUALITY ASSURANCE AND QUALITY CONTROL

The meter(s) should be calibrated before sampling trip following the manufacturer's instructions and the calibrations should be recorded. Additionally, temperature probes will be checked annually against a certified NIST thermometer.

11.0 REFERENCES

Standard Methods.
17th Edition, APHA, 1989.

APPENDIX A
Field Reporting Forms



River and Stream Sampling Field Log
North Dakota Department of Environmental Quality
Division of Water Quality - Watershed Management Program
Telephone: 701-328-6140
Fax: 701-328-6280

Sample #:		Site ID:	Site Description:		Comments:
Dup	Blk	Date: / /	Temperature	DO	
		Time: :	SC	pH	
Sample #:		Site ID:	Site Description:		Comments:
Dup	Blk	Date: / /	Temperature	DO	
		Time: :	SC	pH	
Sample #:		Site ID:	Site Description:		Comments:
Dup	Blk	Date: / /	Temperature	DO	
		Time: :	SC	pH	
Sample #:		Site ID:	Site Description:		Comments:
Dup	Blk	Date: / /	Temperature	DO	
		Time: :	SC	pH	
Sample #:		Site ID:	Site Description:		Comments:
Dup	Blk	Date: / /	Temperature	DO	
		Time: :	SC	pH	
Sample #:		Site ID:	Site Description:		Comments:
Dup	Blk	Date: / /	Temperature	DO	
		Time: :	SC	pH	



Wadable Wetland Field Sampling Form
North Dakota Department of Environmental Quality
Division of Water Quality – Watershed Management Program
Telephone: 701-328-6140
Fax: 701-328-6280

Site ID:	Wetland Name or description:		County:	
Observer:		Date:	Time:	
Aquatic Zone Description:				
Ambient Temp:		Wind Speed:	Wind Direction:	
Cloud Cover	% of 100	Secchi (m)	chlorophyll-a=(y / n)	Phytoplankton=(y / n)
Comments:				
Water Chemistry Taken At		Meters,	Meters,	Meters,

Depth (m) or Location	Temp. (C)	D.O. (mg/L)	pH	Conductivity (umhos/cm)	Comments

Site Drawing:

APPENDIX B
SOP Acknowledgement and Training Form

